

M.E. Schedl's Putnam Lab Open Lab Notebook

Notebook of a Lab Manager for the Putnam Lab

[Notebook](#) [Website](#) [About and CV](#)

browse by [category](#) or [date](#) or [tag](#)

[admin](#)

Methylation Comparison WGBS With Pico Methyl Seq Library Prep Kit First Try

Using the [Zymo Pico Methyl Seq Kit](#) on 3 *Pocillopora* and 3 *Montipora* samples for the purpose of whole genome bisulfate sequencing to compare to RRBS and MBD-BS preparations

Samples for methylation comparison are from the Holobiont Integration Project, and were extracted by Emma Strand or myself, see her notebook posts for the extraction specifications: [20190805](#) and [20180823](#), [20190718](#) and [20190903](#).

Started with 100ng of DNA for each sample, and spiked in 5ng of unmethylated lamda DNA (originally from the RRBS kit) to test bisulfite conversion. Then increaesd the volume to 20µl for the start of bisulfite conversion.

Samples are all Soft homogenization DNA, and if elutions were separated they were the second elution.

Sample	volume DNA (100ng)	Volume spike in	Volume 10mM Tris HCl
--------	--------------------	-----------------	----------------------

Sample	volume DNA (100ng)	Volume spike in	Volume 10mM Tris HCl
1041	1.6	2	16.4
1471	1.8	2	16.18
1637	1.15	2	16.85
1101	3.69	2	14.31
1548	4.33	2	13.67
1628	5.88	2	12.12

Bisulfite Conversion

- Added 130 μ l lightning conversion reagent to each sample in at PCR tube
- Put tubes in the thermocycler Pico bisulfite conversion program

Cleanup

- Added 600 μ l M-binding buffer each to six spin columns
- Added 150 μ l of the BS reaction (all) to each individual tube
- Invert columns to mix
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 100 μ l M-Wash buffer to each column
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 200 μ l L-desulfonation buffer to each column
- Let them sit for 15 minutes
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 8 μ l warmed (56C) DNA elution buffer to each column and let sit for 1 minute in new 1.5mL tubes to collect
- Centrifuged columns at 12,000 rcf for 30 seconds

Amplification with PrepAmp Primers

- Made Priming master mix on ice:
 - $2\mu\text{l}$ 5X PrepAmp buffer * 6.2 = $12.4\mu\text{l}$
 - $1\mu\text{l}$ PrepAmp Primers (40 μM) * 6.2 = $6.2\mu\text{l}$
- Made new PCR tubes with $3\mu\text{l}$ of PrepAmp MM and $7\mu\text{l}$ of bisulfite treated DNA
- Kept those on ice
- Made PrepAmp Mix on ice:
 - $1\mu\text{l}$ 5X PrepAmp buffer * 6.2 = $6.2\mu\text{l}$
 - $3.75\mu\text{l}$ PrepAmp PreMix * 6.2 = $23.25\mu\text{l}$
 - $0.3\mu\text{l}$ PrepAmp polymerase * 6.2 = $1.86\mu\text{l}$
- Set thermocycler program with lid temp restricted to 25 degrees C and place samples inside and run:
 - 98 for 2 minutes
 - 8 degrees for 1 minute
 - 8 degree hold
 - **During hold vortex, spin tubes down, add $5.05\mu\text{l}$ PrepAmp Mix to each tube, vortex, spin down, and place back in thermocycler**
 - 8 degrees for 4 minutes
 - 16 degrees for 1 minute with 3% ramp rate
 - 22 degrees for 1 minute with 3% ramp rate
 - 28 degrees for 1 minute with 3% ramp rate
 - 36 degrees for 1 minute with 3% ramp rate
 - 36.5 degrees for 1 minute with 3% ramp rate
 - 37 degrees for 8 minutes
 - repeat back from the first step one time through again
 - **During hold, vortex, spin tubes down, tried to add $0.3\mu\text{l}$ PrepAmp Polymerase to each tube, it wouldn't really come out of the tip so I ended up adding $.5\mu\text{l}$, vortex, spin down, and place back into thermocycler**

Cleanup with DNA Clean and Concentrator Columns (DCC)

- Made a 1.5mL tube for each sample, added 7:1 ratio DNA binding buffer, so $107.45\mu\text{l}$ of DNA binding buffer

- Put elution buffer in thermomixer 56 degrees
- Added DNA sample (15.35 μ l) to the appropriate 1.5mL tube
- Vortexed, spun down, and added to the column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Transferred columns to 1.5mL tubes
- Added 12 μ l warmed elution buffer to each column directly
- Incubated 1 minute
- Centrifuged 12,000 rcf 30 seconds

First Amplification

- Made 1st Amp master mix:
- 12.5 μ l 2X Library Amp Mix * 4.2 = 52.5 μ l
- 1 μ l Library Amp Primer(10 μ M) * 4.2 = 4.2 μ l
- Added 13.5 μ l MM to new PCR tubes
- Added 11.5 μ l of cleaned and concentrated DNA sample to the appropriate new PCR tube *note here, samples 2 and 6 had a lot more than 11.5 μ l in them, could have been carryover of wash buffer in the column, only 11.5 μ l were used for this step though to keep volumes the same. What was left over I froze in the -20*
- Vortexed, spun down, and placed in thermocycler program 1st Pico Methyl Amp program 6 cycles for the 100ng input

Cleanup with DCC

- Made a 1.5mL tube for each sample, added 7:1 ratio DNA binding buffer, so 175 μ l of DNA binding buffer
- Put elution buffer in thermomixer 56 degrees
- Added DNA sample (25 μ l) to the appropriate 1.5mL tube
- Vortexed, spun down, and added to the column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added 200 μ l M-wash buffer to each column

- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Transferred columns to 1.5mL tubes
- Added 12.5 μ l warmed elution buffer to each column directly
- Incubated 1 minute
- Centrifuged 12,000 rcf 30 seconds

Amplification with Index Primers

- Made new PCR tubes for each sample with 12.5 μ l Library Amp Mix
- **Index primers are the dual index UDI i5 and i7 indexes. That means each sample needs to get primers from 2 tubes, but the total volume of primers needs to be .5 μ l (so .25 μ l each) which is already a lot smaller than what can be easily pipetted. Because of previous issues pipetting .3 μ l, I made separate tubes with 1 μ l of each primer pair and then added .5 μ l of that pre-made pair to each sample**
- Added indexed primers:

sample	index pair
1041	i5 1 and i7 1
1471	i5 2 and i7 2
1637	i5 3 and i7 3
1101	i5 10 and i7 10
1548	i5 11 and i7 11
1628	i5 12 and i7 12



- Added 12 μ l of sample to the appropriate tube (all of the flowthrough from above DCC)
- Vortexed, spun down, and placed in thermocycler program 2nd Pico Methyl Amp program



Cleanup with DCC



- Made a 1.5mL tube for each sample, added 7:1 ratio DNA binding buffer, so 175 μ l of DNA binding buffer
- Put elution buffer in thermomixer 56 degrees
- Added DNA sample (25 μ l) to the appropriate 1.5mL tube
- Vortexed, spun down, and added to the column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added 200 μ l M-wash buffer to each column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Transferred columns to 1.5mL tubes
- Added 12 μ l warmed elution buffer to each column directly
- Incubated 1 minute
- Centrifuged 12,000 rcf 30 seconds

D5000 TapeStation

See the full report [here](#)

B1: 1041 WGBS**C1: 1471 WGBS**

D1: 1637 WGBS**E1: 1101 WGBS**

F1: 1548 WGBS**G1: 1628 WGBS**

Here, 2 and 6 didn't work. Which might be because of the extra volume after the DCC before the first amplification. In this [post](#) which was a test of this kit, one of the samples had extra volume after a cleanup step and then did not have anything in the sample afterwards either. My guess with what is going on with the other ones is that there was too

much input, these traces look similar to 221 NC in that same above post, where I had also started with 100ng input.

Written on September 19, 2019
