

## M.E. Schedl's Putnam Lab Open Lab Notebook

Notebook of a Lab Manager for the Putnam Lab

[Notebook](#) [Website](#) [About and CV](#)

---

browse by [category](#) or [date](#) or [tag](#)

[admin](#)

# RRBS Library Prep of Samples for Methylation Comparison

## Reduced Representation Bisulfite Sequencing Library Prep with the Zymo Research EZ DNA Methylation RRBS Library Prep Test Kit for Methylation Comparison

Samples for methylation comparison are from the Holobiont Integration Project, and were extracted by Emma Strand or myself, see her notebook posts for the extraction specifications: [20190805](#) and [20180823](#), [20190718](#) and [20190903](#).

### Mspl Digestion

- Made PCR tubes with 300ng of DNA, 15ng of Spike-in non-methylated eColi DNA, and water to 35.5 $\mu$ l

Sample	Volume DNA	Volume spike	Volume ultra pure water
1041	4.81	6	24.69
1471	5.47	6	24.03
1637	3.46	6	26.04
1101	11.07	5	18.43
1548	12.99	6	16.51
1628	17.65	6	11.85

- Made master mix for digestion
  - 10X RRBS buffer  $4\mu\text{l} * 6.2 = 24.8\mu\text{l}$
  - MspI  $0.5\mu\text{l} * 6.2 = 3.1\mu\text{l}$
- Added  $4.5\mu\text{l}$  of digestion master mix to each sample and pipetted to mix and spun down
- Placed in thermocycler RRBS digestion program: 4 hours at 37 degrees C and then a 4 degree hold

## Adapter Ligation

- Made master mix for adapter ligation
  - 10X RRBS buffer  $1\mu\text{l} * 6.2 = 6.2\mu\text{l}$
  - rATP  $0.5\mu\text{l} * 6.2 = 3.1\mu\text{l}$
  - MspI  $1\mu\text{l} * 6.2 = 6.2\mu\text{l}$
  - T4 DNA ligase  $1\mu\text{l} * 6.2 = 6.2\mu\text{l}$
  - ultra pure water  $6\mu\text{l} * 6.2 = 37.2\mu\text{l}$
- Pipetted to mix master mix because ligase enzyme
- Added  $9.5\mu\text{l}$  of master mix to each tube
- added  $0.5\mu\text{l}$  RRBS adapters to each tube (adding separately minimizes adapters ligating to themselves)
- Pipetted to mix
- Put in thermocycler program RRBS adapters overnight:

- 21 degrees C 3 hours
- 37 degrees C 1 hour
- 20 degrees C 1 hour
- 37 degrees C 1 hour
- 20 degrees C 1 hour
- 4 degree hold
- The next morning made gap filling master mix
  - Taq DNA Polymerase  $0.5\mu\text{l} * 6.2 = 3.1\mu\text{l}$
  - 5-methylcytosine dNTP Mix  $1.5\mu\text{l} * 6.2 = 9.3\mu\text{l}$
- Added  $2\mu\text{l}$  of gap filling mix to each tube and pipetted to mix
- Placed in thermocycler program 74 degrees C for 30 minutes

### Cleanup with DCC

- Made a 1.5mL tube for each sample, added 7:1 ratio DNA binding buffer, so  $364\mu\text{l}$  of DNA binding buffer
- Put elution buffer in thermomixer 56 degrees
- Added DNA sample ( $52\mu\text{l}$ ) to the appropriate 1.5mL tube
- Vortexed, spun down, and added to the column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added  $200\mu\text{l}$  M-wash buffer to each column
- Centrifuged 12,000 rcf 30 seconds, discarded flowthrough
- Added  $200\mu\text{l}$  M-wash buffer to each column
- Centrifuged 12,000 rcf **1 minute and 30 seconds**, discarded flowthrough
- Transferred columns to 1.5mL tubes
- Added  $20\mu\text{l}$  warmed elution buffer to each column directly
- Incubated 1 minute
- Centrifuged 12,000 rcf 30 seconds

### Bisulfite Conversion

- Transferred each sample to PCR tubes
- Added  $130\mu\text{l}$  lightning conversion reagent to each sample in at PCR tube
- Put tubes in the thermocycler Pico bisulfite conversion program (the same for this kit)

## Cleanup

- Made 6 spin column, one for each sample
- Added 600 $\mu$ l M-binding buffer each to 6 spin columns
- Added 150 $\mu$ l of the BS reaction (all) to each individual tube
- Invert columns to mix
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 100 $\mu$ l M-Wash buffer to each column
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 200 $\mu$ l L-desulfonation buffer to each column
- Let them sit for 15 minutes
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 200 $\mu$ l M-wash buffer to each column
- Centrifuged columns at 12,000 rcf for 30 seconds and discarded flowthrough
- Added 200 $\mu$ l M-wash buffer to each column
- Centrifuged columns at 12,000 rcf for **1 minute and 30 seconds** and discarded flowthrough
- Added 24 $\mu$ l warmed (56C) DNA elution buffer to each column and let sit for 1 minute in new 1.5mL tubes to collect
- Centrifuged columns at 12,000 rcf for 30 seconds

## Index Primer Amplification

- Transferred 24 $\mu$ l of DNA to new PCR tubes
- Added 25 $\mu$ l of Library Amp Master mix to each tube
- Added 0.5 $\mu$ l of each planned primer pairs:

Sample	i5	i7
1041	4	4
1471	5	5
1637	6	6

Sample	i5	i7
1101	13	13
1548	14	14
1628	15	15


- Placed in thermocycler RRBS index amp program for 11 cycles:
  - 94 for 30 sec
  - **94 for 30 sec** (11)
  - **55 for 30 sec** (11)
  - **68 for 1 min** (11)
  - 68 for 5 min
  - 4 hold

### KAPA Pure Bead Cleanup


- Took beads out of fridge beforehand and made fresh 80% EtOH
- Added 50 $\mu$ l (1X) of Kapa Pure beads to each sample
- Performed normal bead cleanup
- Resuspended and eluted samples in 16 $\mu$ l DNA elution buffer

### D5000 TapeStation


### Full Results

**B1: 1041 RRBS****Sample Table**


Well	Conc. [ng/ul]	Sample Description	Alert	Observations
B1	7.78	1041 RRBS		

**C1: 1471 RRBS****Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
C1	3.57	1471 RRBS		


**D1: 1637 RRBS****Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
D1	7.35	1637 RRBS		

**E1: 1101 RRBS****Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
E1	3.14	1101 RRBS		


**F1: 1548 RRBS**



**Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
F1	1.17	1548 RRBS		



**G1: 1628 RRBS****Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
G1	3.73	1628 RRBS		

They have long tails but they all look like they worked!!

*Written on October 24, 2019*